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(一) Mini-M plasmid DNA extraction system (VIOGENE)

Grow 1 to 5ml plasmid-containing bacterial cells in LB medium with appropriate antibiotics overnight (12-16 hours) with vigorous agitation

Pellet the cells by centrifuging for 1-2 minutes. Decant the supernatant and remove all medium residue by pipet.

Add 250 μ l of MX1 Buffer to the pellet, and resuspend the cells completely by vortexing or pipetting.

Add 250 μ l of MX2 Buffer and gently mix (Invert the tube 4 ~ 6 times) to lyse the cells until the lysate becomes Clear. Incubate at room temperature for 1-5 minutes.

Add 350 μ l of MX3 Buffer to neutralize the lysate, then immediately and gently mix the solution. A White precipitate should be formed.

Centrifuge for 5~10 minutes, meanwhile place a mini-M Column onto a Collection Tube.

Transfer the supernatant carefully into the column

Centrifuge for 30-60 seconds. Discard the flow-through.

Wash the column once with 0.5 ml WF Buffer by centrifuging for 30-60 seconds. Discard the flow-through.

Wash the column once with 0.7 ml WS Buffer by centrifuging for 30-60 seconds.
Discard the flow-through.

Centrifuge the column at full speed for another 3 minutes to remove residual ethanol.

Place the column onto a new 1.5-ml centrifuge tube. Add 50µl of Elution Buffer (provided) onto the center of the membrane.


Stand the column for 1-2 minutes, and centrifuge for 1-2 minutes to elute DNA.

Store plasmid DNA at 4°C or -20°C.

(二) Isolation of plasmid DNA (small-scale)(傳統方法)

solution I	solution II	solution III	
50 mM glucose	0.2 N NaOH	5M KoAC	60 ml
25 mM Tris-Cl (pH 8.0)	1% SDS	glacial acetic acid	11.8 ml
10 mM EDTA (pH 8.0)		water	28.5 ml

RNAase(stock solution 10 mg/ml)working solution 50 µg/ml

(三) Genomic DNA extraction miniprep system (VIOGENE)

Pipet up to 200 μ l sample into a 1.5-ml sterile eppendorf tube. When the sample volume is less than 200 μ l, add PBS to make up to 200 μ l. (破細胞的方法因樣品來源而異，植物細胞常用液態氮磨碎，真菌以酵素水解)

Add 20 μ l Proteinase K and 200 μ l EX Buffer into the sample. Mix immediately by vortexing for 20 seconds.

Incubate at 60°C for 20 minutes to lyse the sample. Vortex or invert mix the sample every 3-5 minutes during incubation.

Adjust the incubator to 70°C to incubate for 20 minutes.

Meanwhile, preheat 10 mM Tris-HCl (pH 9.0), ddH₂O, or TE buffer (provided by user) at 70°C (500 μ l /prep) for DNA elution at Step 10.

Add 210 μ l of absolute ethanol or isopropanol to the sample from Step 4 and mix by vortexing.

Place a B/T Genomic DNA Mini Column onto a Collection Tube. Pipette all the mixture (including any precipitate) into the column without touching the rim. Centrifuge at 8000 rpm for 2 minutes. Place the column onto a new Collection Tube.

Wash the column twice with 0.5ml WS buffer by centrifuging at 8000 rpm for 2 minutes. Discard the flow-through after centrifugation.

Centrifuge the column at full speed for another 2 minutes to remove ethanol residue.

Place the column onto a new 1.5 ml tube (provided by user). Elute DNA with 200 μ l of the preheated elution solution from step 5


Stand the column for 1-5 minutes, and centrifuge for 1-2 minutes to elute DNA

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DNA staining solution (ethidium bromide (10 mg/ml), SYBR Gold)

Electrophoresis buffer (usually 1x TAE or 0.5x TBE)

6X Gel-loading buffer (0.25% (w/v) bromophenol blue, 40% (w/v) sucrose in H₂O)

Nucleic Acids and Oligonucleotides

DNA samples

DNA size standards

Samples of DNAs of known size are typically generated by restriction enzyme digestion of a plasmid or bacteriophage DNA of known sequence. Alternatively, they are produced by ligating a monomer DNA fragment of known size into a ladder of polymeric forms.

Agarose Concentration in Gel (% [w/v])	Range of Separation of Linear DNA Molecules (kb)
0.3	5-60
0.6	1-20
0.7	0.8-10
0.9	0.5-7
1.2	0.4-6
1.5	0-2-3
2.0	0.1-2

Agarose gels are cast by melting the agarose in the presence of the desired buffer until a clear, transparent solution is achieved. The melted solution is then poured into a mold and allowed to harden. Upon hardening, the agarose forms a matrix, the density of which is determined by the concentration of the agarose.

Ethidium Bromide

Add 1 g of ethidium bromide to 100 ml of H₂O. Stir on a magnetic stirrer for several hours to ensure that the dye has dissolved. Wrap the container in aluminum foil or transfer the 10 mg/ml solution to a dark bottle and store at room temperature.

TAE

Prepare a 50X stock solution in 1 liter of H₂O:

242 g of Tris base

57.1 ml of glacial acetic acid

100 ml of 0.5 M EDTA (pH 8.0)

The 1x working solution is 40 mM Tris-acetate/1 mM EDTA.

TBE

Prepare a 5X stock solution in 1 liter of H₂O:

54 g of Tris base

27.5 g of boric acid

20 ml of 0.5 M EDTA (pH 8.0)

The 0.5x working solution is 45 mM Tris-borate/1 mM EDTA.

EDTA

To prepare EDTA at 0.5 M (pH 8.0): Add 186.1 g of disodium EDTA•2H₂O to 800 ml of H₂O. Stir vigorously on a magnetic stirrer. Adjust the pH to 8.0 with NaOH (approx. 20 g of NaOH pellets). Dispense into aliquots and sterilize by autoclaving. The disodium salt of EDTA will not go into solution until the pH of the solution is adjusted to approx. 8.0 by the addition of NaOH.

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

High GC content target need to add DMSO

PCR PROGRAM

動作	Stage	Step	溫度	時間	Cycle
前置工作	1	1	94°C	10min	1
Denature	2	1	94°C	1 min	30
Annealing	2	2	64°C	35sec	
Extension	2	3	72°C	45sec	
Extra extension	3	1	72°C	7min	1

(七) PCR-M Clean-up system kit (PCR 產物純化 KIT)

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